# OW-SERVER

[www.AgileBio.com](http://www.AgileBio.com)

---

## 1 Installation d'un OW-Server

### 1.1 Branchement du OW-Server-Enet-2

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
\*


- WEP\_40: Wired Equivalent Privacy, requires 5 ASCII characters or 10 hexadecimal characters (0-9 or A-F)
- WEP\_104: Wired Equivalent Privacy, requires 13 ASCII characters or 26 hexadecimal characters (0-9 or A-F).
- WPA\_Key: WiFi Protected Access, also called WPA-Personal. Requires 64 hexadecimal characters (0-9 or A-F).
- WPA\_Auto\_Key: Device will connect to access point using WPA2-Personal if it is supported, otherwise it will use WPA-Personal. Requires 64 hexadecimal characters (0-9 or A-F).
- WPA\_Passphrase: WiFi Protected Access, also called WPA-Personal. Requires between 8 and 63 ASCII characters.
- WPA\_Auto\_Passphrase: Device will connect to access point using WPA2-Personal if it is supported, otherwise it will use WPA-Personal. Requires between 8 and 63 ASCII characters.
- WPA2\_Key: WiFi Protected Access II, also called WPA2-Personal. Requires 64 hexadecimal characters (0-9 or A-F).
- WPA2\_Passphrase: WiFi Protected Access II, also called WPA2-Personal. Requires between 8 and 63 ASCII characters.

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
WLAN Type:  AdHoc  
 Infrastructure  
SSID Name: EDSWIFI  
Security Type: Open  
Security Key:  
Save

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
System Configuration ▾

- Network
- WiFi
- POST Client
- 1-Wire Interface

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




A screenshot of a network configuration form. The "Enable DHCP" checkbox is checked. The form contains the following fields:

IP Address:	169.254.1.1
Gateway:	169.254.1.1
Subnet Mask:	255.255.0.0
Primary DNS:	169.254.1.1
Secondary DNS:	169.254.1.1
HTTP Port:	80
Save	

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
The screenshot shows the OW Server v2 Enet web interface. At the top left is the logo for EDS EMBEDDED DATA SYSTEMS. At the top right, it says "OW Server v2 Enet" and "OWServer\_v2-Enet". Below the logo is a navigation menu with "Home", "Advanced", "System Configuration", and "Contact us". The main content area displays system status: "Devices Connected: 1", "Loop Time: 1.037 Sec", "Poll Count: 10141300", and "Supply Voltage: 5.13". It also shows data for three channels: Channel 1 (Devices: 1, Errors: 0, Voltage: 4.90), Channel 2 (Devices: 0, Errors: 0, Voltage: 4.90), and Channel 3 (Devices: 0, Errors: 0, Voltage: 4.89). A table below shows a "Description" of a "Programmable resolution thermometer" with a "ROM ID" of "70000003E6B0B228", "Device" of "DS18B20", "Channel" of "1", "Health\*" of "7", and "Value" of "75.9 °F". At the bottom, there are dropdown menus for "Duration" (15 Minutes), "Channel" (All), and "Value" (None), along with "Begin" and "Download" buttons.

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The screenshot shows the "POST Client Configuration" page. It starts with the title "POST Client Configuration" and a subtitle "Use this page to configure the automatic sending of data to an HTTP server using POST." Below this, there are two sections: "Enable Client:" with a checked checkbox and "URL:" with the text "http:// LABCOLLECTOR URL /UserFiles/lclogger\_autoload\_ow.php". The "Period:" is set to "60" seconds. The "Enable Proxy:" section has an unchecked checkbox and "Proxy URL:" set to "http://myproxy.com:8080". A "Save" button is at the bottom of the configuration area. Below the configuration area, there is a "Notes:" section. The first note says "The latest device data can be automatically sent to a HTTP Post server using the POST Client feature. Data is sent in XML format." The second note says "URL - Where to post the data, as typed into a browser. Examples:" followed by three examples: "http://www.embeddeddatasystems.com/upload.htm", "http://www.embeddeddatasystems.com:12391/files/myfile.htm", and "http://192.168.1.23/upload.htm". The third note says "Period - How often, in seconds, to send the data." The final note says "HTTP Post can be sent through a proxy server by enabling it and entering the proxy server location in the Proxy URL field using the same rules as above."