LSM Lab Service Manager



LSM Add-on brings a complete LIMS/LIS to the LabCollector platform, fully and easily configurable to fit all types of analytical and other service labs.












Core labs

Simplifying the process from ordering to reporting























♦ Streamlined service delivery **♦** Fully configurable

Comprehensive job and workload overview

Sample/batch processing with barcode support

100% Web and mobile/tablet compatible

Integrated billing system and reporting tool

Requester address book & shipment manifests

Integration with ERP/CRM (SAP, Dynamics...) via API

Easy test and service configuration

Report and invoice template setup and options

Flexible test pricing

Customizable data fields and forms

Dynamic data import and export capabilities

HL7, ASTM, CSV, XMC... integration

≝Compliancy

ISO 17015 compliant: full audit trail / result validation

Electronic signatures and authentication

Document/data control and versioning

User activity logs and monitoring

Data trending and Control charts

Security

Credentials with possibility of LDAP/AD and SSO

Multiple user permission levels and groups

Data Encryption and backup

Regular security audits and pentests

Secure API

Ready to use remote web portal

Requesters, clients or partners can submit requisitions, view order status, download PDF Reports with results, Quotes and Invoices

Ready-to-use 100 % web extension of LSM Add-on

Easy to configure web portal, rebrand and embed in your company website


Patient Portal mode with kit activation



Automation & Instrument Integration

Connect any instrument and build automation workflows to improve efficiency and data integrity.

Any type of instrument of any brand: Analyzers, Diagnostics and Chromatography Systems ... Easy HL7/ASTM/XMC/CSV, data mapping.



AgileBio LLC San Diego, USA Phone: +1 347 368-1315 sales.us@agilebio.com

AgileBio SARL Phone: +33 1 41 79 15 85 sales@agilebio.com

