



DATA LOGGER

LabCollector LIMS Add-on

Laboratory Storage Temperature Tracking



Laboratory Equipment



Data logger and temperature sensors



Data logger add-on for LabCollector

SMS Alerts!



GSM box for SMS alerts


- Network based application for day-to-day control of storage equipments temperature, humidity, CO₂ and other data centralized in LabCollector LIMS.
- Temperature log visible for each sample in the database. Retrieve data for a defined period (quality control). Historical charts available.
- Email and/or SMS alerts when values are not correct.
- Low temperature sensors (-80°C and -200°C).

Contact USA
 sales.us@agilebio.com
 Phone: +1 (347) 368-1315
 Fax: +1 (800) 453-9128

LabCollector
 www.LabCollector.com

Contact France and Worldwide
 sales@agilebio.com
 Phone: +33 (0)1 41 79 15 85
 Fax: +33 (0)1 72 70 40 22

Easy plug and play metrology system for your lab



▶ Data logger to connect your sensors: temperature, humidity, CO₂... (Ethernet/Wifi).



▶ Analogic extensions for RTD sensors (-200°C).



▶ GSM/SMS box for sending alerts.



▶ RJ45 extension.



▶ Low temperature RTD sensors (-200°C/100°C).



▶ Digital temperature sensors (-55°C/150°C).

Ask your FREE* testing kit!



*To allow testing this add-on in optimal conditions, we can send you a 30 days free testing kit (1 Ethernet logger + 1 digital temperature sensor)